

Product Datasheet

Human SERPINB5 protein (orb1463669)

Description	Recombinant Human SERPINB5 protein
Reactivity	Human
Tested Applications	ELISA, WB
Preservatives	Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose and 0.06% proclin, pH7.4
Storage	-20°C for 12 months as lyophilized; 2-8°C for 1 month under sterile conditions after reconstitution
Tag	N-terminal His-IF2DI Tag
Note	For research use only
Application notes	Centrifuge the vial at 10,000 rpm for 1 minute, reconstitute at 200 µg/ml in sterile distilled water by gentle pipetting 2-3 times, don't vortex
Formula	Lyophilized from a 0.2 µm filtered solution in 10 mM Hepes, 500 mM NaCl with 5% trehalose and 0.06% proclin, pH _{7.4}
Purity	Greater than 90% by SDS-PAGE gel analyses
Source	E.Coli
MW	50.4 kDa
Uniprot ID	P36952
Expiration Date	6 months from date of receipt.