

Product Datasheet

Mouse IL-36 gamma Protein (orb1463839)

Description	The Mouse IL-36 gamma yeast-derived recombinant protein is not tagged and is naturally endotoxin-free...
Target	IL-36 gamma
Form/Appearance	Lyophilized
Storage	-20°C
Note	For research use only
Protein Sequence	GRETPDFGEV FDLQQVWIF RNQALVTVPR SHRVTPVSVT ILPCKYPESL EQDKGIAIYL GIQNPDKCLF CKEVNGHPTL LLKEEKILD YHHPEPMKPF LFYHTRTGGT STFESVAFPG HYIASSKTGN PIFLTSKKGE YYNINFNLDI K (151)
Purity	98%
Source	Yeast
MW	17.2 kDa
Entrez	215257
Expiration Date	6 months from date of receipt.