



Product Datasheet

p70 S6 Kinase (Ab-411) Antibody (orb14646)


| | |
|----------------------------|---|
| Description | p70 S6 Kinase (Ab-411) Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IF, IHC, WB |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | P23443 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of p70 S6 Kinase(A...



Immunohistochemical analysis of formalin...



Confocal immunofluorescent analysis of M...