


Product Datasheet

Human Chitinase 3-like 1/YKL-40 ELISA Kit (orb146700)

Description	Human Chitinase 3-like 1/YKL-40 ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Y22-T383
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:50, concentration in serum and plasma is 12.72-112.2 ng/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and urine.
MW	42625 MW
Uniprot ID	P36222
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human Chitinase 3-like 1/YKL-40 PicoKine...