



## Product Datasheet

### PKCδ (Ab-645) Antibody (orb14671)


<b>Description</b>	PKCδ (Ab-645) Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, IHC, WB
<b>Storage</b>	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<a href="#">Q05655</a>
<b>Expiration Date</b>	12 months from date of receipt.



Confocal immunofluorescent analysis of H...



Western blot analysis of PKCδ (Ab-645) an...



Immunohistochemical analysis of formalin...