



Product Datasheet

SHP (Ab-536) Antibody (orb14690)


Description	SHP (Ab-536) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P29350
Expiration Date	12 months from date of receipt.



Western blot analysis of SHP-1(Ab-536) a...



Confocal immunofluorescent analysis of M...



Immunohistochemical analysis of formalin...