



Product Datasheet

HDAC4/HDAC5/HDAC9 (Ab-246/259/220) Antibody (orb14717)


Description	HDAC4/HDAC5/HDAC9 (Ab-246/259/220) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P56524
Expiration Date	12 months from date of receipt.



Western blot analysis of HDAC4/HDAC5/HDAC9



Confocal immunofluorescent analysis of HDAC4/HDAC5/HDAC9



Immunohistochemical analysis of formalin-fixed tissue