

Product Datasheet

Human CNTF Protein (orb1471705)

Description	Recombinant Human Ciliary Neurotrophic Factor is produced by our E.coli expression system and the ta...
Reactivity	Human
Conjugation	Unconjugated
Target	CNTF
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM Tris-HCl, 100mM NaCl, pH 8.0.
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of 20mM Tris-HCl, 100mM NaCl, pH 8.0.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Note	For research use only
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	E. Coli
MW	22.93 KDa
Uniprot ID	P26441
Expiration Date	6 months from date of receipt.