

Product Datasheet

Human IL-15RAandIL-15 (C-Fc) Protein (orb1471717)

Description	Recombinant Human Interleukin-15 Receptor Alpha and Interleukin-15 Fusion Protein is produced by our...
Reactivity	Human
Conjugation	Unconjugated
Target	IL-15RA&IL-15
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, 5% Trehalose, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of PBS, 5% Trehalose, pH 7.4.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Tag	C-Fc Tag
Note	For research use only
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	Mammalian
MW	46.9 KDa
Uniprot ID	Q13261 , P40933
Expiration Date	6 months from date of receipt.