

## Product Datasheet

### Human IL-22RA2 (C-Fc) Protein (orb1471718)

<b>Description</b>	Recombinant Human Interleukin-22 Receptor Subunit Alpha-2/Interleukin-22 Binding Protein is produced...
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	IL-22RA2
<b>Preservatives</b>	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
<b>Form/Appearance</b>	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.4.
<b>Storage</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Tag</b>	C-Fc Tag
<b>Note</b>	For research use only
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE.
<b>Source</b>	Mammalian
<b>MW</b>	51.8 KDa
<b>Uniprot ID</b>	<a href="#">Q969J5</a>
<b>Expiration Date</b>	6 months from date of receipt.