

Product Datasheet

Human/Mouse FGF-8b Protein (orb1471729)

Description	Recombinant Human/Mouse Fibroblast Growth Factor 8B is produced by our E.coli expression system and ...
Reactivity	Human, Mouse
Conjugation	Unconjugated
Target	FGF-8b
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Note	For research use only
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	E. Coli
MW	22.5 KDa
Uniprot ID	P37237, P55075
Expiration Date	6 months from date of receipt.