

Product Datasheet

Human pro-BDNF Protein (orb1471731)

Description	Recombinant Human Pro-Brain-Derived Neurotrophic Factor is produced by our E.coli expression system ...
Reactivity	Human
Conjugation	Unconjugated
Target	pro-BDNF
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, pH 8.0.
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of PBS, pH 8.0.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Note	For research use only
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	E. Coli
MW	25.6 KDa
Uniprot ID	P23560
Expiration Date	6 months from date of receipt.