

Product Datasheet

Human CCL8 (C-6His) Protein (orb1471733)

Description	Recombinant Human C-C Motif Chemokine 8 is produced by our Mammalian expression system and the target...
Reactivity	Human
Conjugation	Unconjugated
Target	CCL8
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, pH 7.4.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Tag	C-6×His Tag
Note	For research use only
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	Mammalian
MW	9.95 KDa
Uniprot ID	AAI26243.1
Expiration Date	6 months from date of receipt.