

Product Datasheet

Human FGF-4 (136AA) Protein (orb1471743)

Description	Recombinant Human Fibroblast Growth Factor 4 is produced by our E.coli expression system and the tar...
Reactivity	Human
Conjugation	Unconjugated
Target	FGF-4
Preservatives	Lyophilized from a 0.2 µm filtered solution of PBS, 5% Trehalose, pH 7.4.
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of PBS, 5% Trehalose, pH 7.4.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Note	For research use only
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	E. Coli
MW	15.1 KDa
Uniprot ID	P08620
Expiration Date	6 months from date of receipt.