

## Product Datasheet

### Human IL-29 (C-10His) Protein (orb1471744)

|                        |   |
|------------------------|---|
| <b>Description</b>     | Recombinant Human Interferon Lambda-1 is produced by our Mammalian expression system and the target ...   |
| <b>Reactivity</b>      | Human   |
| <b>Conjugation</b>     | Unconjugated  |
| <b>Target</b>          | IL-29   |
| <b>Preservatives</b>   | Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.   |
| <b>Form/Appearance</b> | Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.   |
| <b>Storage</b>         | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| <b>Tag</b>             | C-10×His Tag  |
| <b>Note</b>            | For research use only   |
| <b>Purity</b>          | The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.  |
| <b>Source</b>          | Mammalian   |
| <b>MW</b>              | 21.4 KDa  |
| <b>Uniprot ID</b>      | <a href="#">Q8IU54</a>  |
| <b>Expiration Date</b> | 6 months from date of receipt.  |