

## Product Datasheet

# Human IL-2 Superkine (Mutant, C-6His) Protein (orb1471752)

<b>Description</b>	Recombinant Human Interleukin-2 Superkine is produced by our Mammalian expression system and the tar...
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	IL-2
<b>Preservatives</b>	Lyophilized from a 0.2 µm filtered solution of 4mM HCl.
<b>Form/Appearance</b>	Lyophilized from a 0.2 µm filtered solution of 4mM HCl.
<b>Storage</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Tag</b>	C-6×His Tag
<b>Note</b>	For research use only
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE.
<b>Source</b>	Mammalian
<b>MW</b>	16.5 KDa
<b>Uniprot ID</b>	<b>P60568</b>
<b>Expiration Date</b>	6 months from date of receipt.