

Product Datasheet

Human LR3-IGF-1 Protein (orb1471762)

Description	Recombinant Human LR3 Insulin-Like Growth Factor-I is produced by our E.coli expression system and t...
Reactivity	Human
Conjugation	Unconjugated
Target	LR3-IGF-1
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM NaAc-HAc, 4% Mannitol, pH 4.5.
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of 20mM NaAc-HAc, 4% Mannitol, pH 4.5.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Note	For research use only
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	E. Coli
MW	9.1 KDa
Uniprot ID	P05019
Expiration Date	6 months from date of receipt.