

Product Datasheet

Human/Murine/Rat BDNF Protein (orb1471779)

Description	Recombinant Human Brain-Derived Neurotrophic Factor is produced by our E.coli expression system and ...
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Target	BDNF
Preservatives	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 250mM NaCl, pH 7.2.
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 250mM NaCl, pH 7.2.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Note	For research use only
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	E. Coli
MW	13 KDa
Uniprot ID	P23560
Expiration Date	6 months from date of receipt.