

Product Datasheet

Human/Mouse/Rat BMP-2 Protein (orb1471785)

Description	Recombinant Human/Mouse/Rat Bone Morphogenetic Protein 2 is produced by our E.coli expression system...
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Target	BMP-2
Preservatives	Lyophilized from a 0.2 µm filtered solution of 10mM HAc-NH ₄ Ac, 4% D-Mannitol, pH 4.0.
Form/Appearance	Lyophilized from a 0.2 µm filtered solution of 10mM HAc-NH ₄ Ac, 4% D-Mannitol, pH 4.0.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Note	For research use only
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	E. Coli
MW	13.3 KDa
Uniprot ID	P12643
Expiration Date	6 months from date of receipt.