

## Product Datasheet

### Human G-CSF Protein (orb1471786)

<b>Description</b>	Recombinant Human Granulocyte Colony-Stimulating Factor is produced by our E.coli expression system ...
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Target</b>	G-CSF
<b>Preservatives</b>	Lyophilized from a 0.2 µm filtered solution of 10mM HAc-NaAc, 150mM NaCl, 0.004% Tween 80, 5% Mannitol, pH 4.0.
<b>Form/Appearance</b>	Lyophilized from a 0.2 µm filtered solution of 10mM HAc-NaAc, 150mM NaCl, 0.004% Tween 80, 5% Mannitol, pH 4.0.
<b>Storage</b>	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
<b>Note</b>	For research use only
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE.
<b>Source</b>	E. Coli
<b>MW</b>	18.8 KDa
<b>Uniprot ID</b>	<b>P09919</b>
<b>Expiration Date</b>	6 months from date of receipt.