

Product Datasheet

HIV-1 p24 antibody (orb1471808)

| | |
|--------------------------|--|
| Description | Mouse monoclonal antibody to HIV-1 p24 |
| Reactivity | Virus |
| Preservatives | *HIV-1 p24 coating antibody contains PBS, pH 7.4. * HIV-1 p24 conjugating antibody contains PBS, pH 7.4. |
| Form/Appearance | 2 vials of sterile filtered clear colorless solution. |
| Storage | Stability: HIV-1 p24 antibody although stable at 4°C for 1 week, should be stored below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles |
| Note | For research use only |
| Application notes | Applications: Lateral flow immunoassay |
| Clonality | Monoclonal |
| Purity | Greater than 95%. |
| Expiration Date | 12 months from date of receipt. |