

Product Datasheet

HBe VLP antibody (orb1471809)

Description	Mouse monoclonal antibody to HBe VLP
Reactivity	Virus
Conjugation	Unconjugated
Preservatives	Phosphate buffered saline pH 7.4.
Form/Appearance	sterile filtered clear colorless solution.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Applications: ELISA, Lateral flow immunoassay
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.