

Product Datasheet

HB-Pre-sAg antibody (orb1471810)

Description	Mouse monoclonal antibody to HB-Pre-sAg
Reactivity	Virus
Conjugation	Unconjugated
Immunogen	HBsAg purified from human sera pool
Preservatives	1mg/ml in PBS.
Form/Appearance	Sterile Filtered solution.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Monoclonal Antibodies, Anti-Viral
Isotype	Mouse IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.