

Product Datasheet

Anhui 2013 Protein (orb1471817)

Description	Anhui H7N9 produced in Sf9 Baculovirus cells is a single, glycosylated polypeptide chain containing ...
Preservatives	The Anhui H7N9 solution (1mg/ml) contains 10% Glycerol and Phosphate-Buffered Saline (pH 7.4).
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Viral Antigens- Influenza Virus
Protein Sequence	ADLDKICLGH HAVSNGTKVN TLTERGVEVV NATETVERTN IPRICSKGKR TVDLGQCGLL GTITGPPQCD QFLEFSADLI IERREGSDVC YPGKVFNEEA LRQILRESGG IDKEAMGFTY SGIRTNGATS ACRRSGSSFY AEMKWLLSNT DNAAFPQMTK SYKNTRKSPA LIVWGIHHSV STAEQTKLYG SGNKLVTVGS SNYQQSFVPS PGARPQVNGL SGRIDFHWM LNPNDVTFS FNGAFIAPDR ASFLRGKSMG IQSGVQVDAN CEGDCYHSGG TIISNLPFQN IDSRAVGKCP RYVKQRSLLL ATGMKNVPEI PKGRHHHHHHH
Purity	Greater than 95.0% as determined by SDS-PAGE.
Source	Sf9, Baculovirus cells
Expiration Date	12 months from date of receipt.