

Product Datasheet

VWA2 Protein (orb1471834)

Description	VWA2 Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing...
Preservatives	The VWA2 protein solution contains 20mM Tris-HCl, pH 8.0, 0.8M Urea & 50% glycerol.
Form/Appearance	Sterile Filtered colorless liquid formulation.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. Please avoid freeze thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Purity	Greater than 95% as determined by SDS-PAGE.
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.