
Product Datasheet

CD209 Protein (orb1471871)

Description

CD209 Human Recombinant produced in Sf9

polypeptide chai...

Preservatives

The CD209 solution (0.25mg/ml) contains 10% glycerol and Phosphate-Buffered Saline (pH 7.4).

Form/Appearance

Sterile filtered colorless solution.

Storage

Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles

Note

For research use only

Application notes

Recombinant & Natural Proteins

Protein Sequence

ADPVSKVPSS ISQEQSRQDA IYQNLTQLKA
AVGELSEKSK LQEIQELTQ LKAAVGELPE
KSKLQEIQE LTRLKAAVGE LPEKSKLQEI
YQELTWLAAA VGELPEKSKM QEIQELTRL
KAAVGELPEK SKQEQIQEL TRLKAAVGEL
PEKSKQEQEY QELTRLKAAV GELPEKSKQQ
EIQELTQLK AAVERLCHPC PWEWTFQGN
CYFMSNSQRN WHDSITACKE VGAQLVVIKS
AEEQNFLQLQ SSRSNRFTWM GLSDLNQEGT
WQWVDGSPLL PSFKQYWNRG EPNNVGEEDC
AEFSGNGWND DKCNLAKFWI CKKSAASCSR
DEEQFLSPAP ATPNPPPAVE PKSCDKTHTC
PPCPAPELLG GPSVFLFPPK PKDTLMISRT
PEVTCVVVDV SHEDPEVKFN WYVDGVEVHN
AKTKPREEQY NSTYRVVSVL TVLHQDWLNG
KEYKCKVSNK ALPAPIEKTI SKAKGQPREP
QVYTLPPSRD ELTKNQVSLT CLVKGFYPSD
IAVEWESNGQ PENNYKTPP VLDSGDSFFL
YSKLTVDKSR WQQGNVFSCS VMHEALHNNHY
TQKSLSLSPG KHHHHHHH

Purity

Greater than 90.0% as determined by SDS-PAGE.

Source

Sf9, Baculovirus cells

Expiration Date

12 months from date of receipt.