

## Product Datasheet

### KLK3 Native Protein (orb1471895)

<b>Description</b>	Human Kallikrein-3 produced in Human seminal fluid having a molecular mass of approximately 30kD.
<b>Preservatives</b>	The protein solution (0.2 µm filtered) is in 0.09% NaN <sub>3</sub> , 0.05M phosphate buffer, 150mM NaCl, pH 7.5.
<b>Form/Appearance</b>	Sterile Filtered clear solution.
<b>Storage</b>	Stability: Human Kallikrein-3 should be stored at 2-8°C. Do not freeze!
<b>Note</b>	For research use only
<b>Application notes</b>	Human Virus Test: Starting material donor has been tested and certified negative for antibodies to HIV-1, HIV-2, HCV, HBSAG, Syphilis and HIV/HBV/HCV (PCR)
<b>Purity</b>	Greater than 96.0%.
<b>Source</b>	Human seminal fluid
<b>Expiration Date</b>	12 months from date of receipt.