

Product Datasheet

PLGF Protein (orb1471913)

Description	PLGF Human Recombinant produced in HEK293 cells is a single, glycosylated polypeptide chain (a.a 21-...
Preservatives	PLGF protein (0.25mg/ml) contains 10% glycerol and Phosphate-Buffered Saline (pH 7.4).
Form/Appearance	Sterile filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks.Store, frozen at -20°C for longer periods of time.For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Cytokines And Growth Factors
Protein Sequence	DGSMVPPQQ WALSAAGNGSS EVEVVPFQEV WGRSYCRALRLVDVVSEYPSVEVHMFSPS CVSLLRCTGC CGDENLHCVP VETANVTMQL LKIRSGDRPS YVELTFSQHV RCECRPLREK MKPERRRPPKG RGKRRREKQR PTDCHLCGDA VPRRHHHHHHH
Purity	Greater than 95.0% as determined by SDS-PAGE.
Source	HEK293 cells
Expiration Date	12 months from date of receipt.