

Product Datasheet

k9IL17A Protein (orb1471919)

Description	IL17A Canine Recombinant produced in HEK293 cells is a single, glycosylated polypeptide chain (29-15...
Preservatives	IL17A protein (0.5mg/ml) contains 10% glycerol and Phosphate-Buffered Saline (pH 7.4).
Form/Appearance	Sterile filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks.Store, frozen at -20°C for longer periods of time.For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA).Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Cytokines And Growth Factors
Protein Sequence	FPQNPGCRNT EDKNFPQHVK VNLNILNRNT NSRRPSDYIN RSTSPWNLHR NEDPERYPSV IWEAKCRHLG CVNNEGNINY HMNSVPIQQE ILVLRRESQH CPHSFRLEKM LVAVGCTCVT PIVRHVAHHH HHH
Purity	Greater than 90.0% as determined by SDS-PAGE.
Source	HEK293 cells
Expiration Date	12 months from date of receipt.