

Product Datasheet

k9IL5 Protein (orb1471920)

Description	IL5 Canine produced in Sf9 Baculovirus cells is a single, glycosylated polypeptide chain containing ...
Preservatives	The IL5 solution (0.25mg/ml) contains 10% Glycerol and Phosphate-Buffered Saline (pH 7.4).
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Cytokines And Growth Factors
Protein Sequence	VENPMNRLVA ETLTLLSTHR TWLIGDGNLM IPTPENKNHQ LCIKEVFQGI DTLKNQTAHG EAVDKLFQNL SLIKEHIERQ KKRCAGERWR VTKFLDYLVQ FLGVINTEWT PESHHHHHH
Purity	Greater than 90.0% as determined by SDS-PAGE.
Source	Sf9, Baculovirus cells
Expiration Date	12 months from date of receipt.