

Product Datasheet

His Tag antibody (orb1472074)

Description	Mouse monoclonal antibody to His Tag
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Synthesized peptide of His Tag (AA: HHHHHHHH-KLH).
Preservatives	Purified antibody in PBS with 0.05% sodium azide
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Dilution Range	ELISA: 1/10000
Expiration Date	12 months from date of receipt.