



Product Datasheet

Synaptophysin Like Protein 1 Rabbit mAb Antibody (orb1472567)

Description	Synaptophysin Like Protein 1 Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q16563
Expiration Date	12 months from date of receipt.



Immunohistochemistry of SYPL1 in paraffi...



Western blot detection of K562 cell lysa...