



Product Datasheet

NgR3 Rabbit mAb Antibody (orb1472609)

Description	NgR3 Rabbit mAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Preservatives	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% rAlbumin
Concentration	0.3 mg/ml
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
Purity	Affinity Purified
MW	Calculated MW:49 kDa, Observed MW:49 kDa
Uniprot ID	Q86UN2
Expiration Date	12 months from date of receipt.



Western blot detection of NgR3 in A549, ...



Immunocytochemistry of NgR3 (green) in H...