



Product Datasheet

CDK6 (Ab3) Antibody (orb14728)

Description	CDK6 (Ab3) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q00534
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Confocal immunofluorescent analysis of H...