

## Product Datasheet

### Sivopixant (orb1473203)

<b>Description</b>	Sivopixant (S-600918) is a potent and selective P2X3 receptor antagonist (P2X3 IC50=4.2 nM; P2X2/3 I...
<b>Note</b>	For research use only
<b>Formula</b>	C <sub>25</sub> H <sub>22</sub> CLN <sub>5</sub> O <sub>5</sub>
<b>MW</b>	507.93
<b>CAS Number</b>	[2414285-40-6]
<b>SMILES</b>	<chem>C(C1C=CC(CI)=CC=1)N1C(N(C[C@H](C)C(=O)O)C(=O)N=C1NC1C=CC(OC2N=CC=CC=2)=CC=1)=O</chem>
<b>Expiration Date</b>	12 months from date of receipt.