


---

## Product Datasheet

**COVID-19 S Protein RBD (Omicron, BA.2)  
(orb1473268)**

**Description**

Recombinant COVID-19 S Protein RBD (Omicron,


**Reactivity**

Virus

**Endotoxins**

Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.

**Tested Applications**

ELISA, MS, SDS-PAGE, WB

**Preservatives**

Liquid in a 0.2 μm filtered solution of PBS, pH 7.4.

**Form/Appearance**

Liquid in a 0.2 μm filtered solution of PBS, pH 7.4.

**Storage**

Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

**Note**

For research use only

**Application notes**

Please minimize freeze-thaw cycles.

**Purity**

Greater than 95% as determined by reducing SDS-PAGE.


**Source**

Human cells


**Uniprot ID**
[P0DTC2](#)
**Entrez**
[43740568](#)
**Expiration Date**

6 months from date of receipt.

Greater than 95% as determined by reduci...



Loaded Human ACE-2 on Protein A Biosenso...



Loaded Anti-COVID-19 S1 mAb on Protein A...