

## Product Datasheet

# COVID-19 S Protein RBD (Omicron, BA.2) (orb1473269)

<b>Description</b>	Recombinant COVID-19 S Protein RBD (Omicron, BA.2)
<b>Reactivity</b>	Virus
<b>Endotoxins</b>	Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.
<b>Tested Applications</b>	ELISA, MS, SDS-PAGE, WB
<b>Preservatives</b>	Liquid in a 0.2 μm filtered solution of PBS, pH 7.4.
<b>Form/Appearance</b>	Liquid in a 0.2 μm filtered solution of PBS, pH 7.4.
<b>Storage</b>	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Please minimize freeze-thaw cycles.
<b>Purity</b>	Greater than 95% as determined by reducing SDS-PAGE.
<b>Source</b>	Human cells
<b>Uniprot ID</b>	<a href="#">P0DTC2</a>
<b>Entrez</b>	<a href="#">43740568</a>
<b>Expiration Date</b>	6 months from date of receipt.