


Product Datasheet

COVID-19 Nucleocapsid Protein (Omicron, B.1.1.529) (orb1473270)

Description	Recombinant COVID-19 Nucleocapsid Protein (Omicron, B.1.1.529)
Reactivity	Virus
Endotoxins	Endotoxin: Less than 0.1 ng/μg (1 IEU/μg) as determined by LAL test.
Tested Applications	ELISA, MS, SDS-PAGE, WB
Preservatives	Liquid in a 0.2 μm filtered solution of 20 mM Tris-HCl, 300 mM NaCl, 1 mM EDTA, pH 8.0.
Form/Appearance	Liquid in a 0.2 μm filtered solution of 20 mM Tris-HCl, 300 mM NaCl, 1 mM EDTA, pH 8.0.
Storage	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Please minimize freeze-thaw cycles.
Purity	Greater than 95% as determined by reducing SDS-PAGE.
Source	E.coli
Uniprot ID	P0DTC9
Entrez	43740575
Expiration Date	6 months from date of receipt.



Greater than 95% as determined by reduci...