


Product Datasheet

Rat Luteinizing Hormone ELISA Kit (orb1473274)

Description	Rat Luteinizing Hormone ELISA Kit
Reactivity	Rat
Range	0.468 ng/ml - 30 ng/ml
Tested Applications	EIA
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma or Other Biological Fluids
Sample Types	Serum, Plasma or Other Biological Fluids
Uniprot ID	P01230
Sensitivity	0.14 ng/ml
Expiration Date	Please enquire.



Representative
Standard
Curve