


## Product Datasheet

### Mouse FGF15 ELISA Kit (orb1473276)

<b>Description</b>	Mouse FGF15 ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	31.2 pg/ml - 2000 pg/ml
<b>Tested Applications</b>	EIA
<b>Target</b>	FGF15
<b>Storage</b>	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Tissue Homogenates, Cell Lysates, Cell Culture Supernates or Other Biological Fluids
<b>Sample Types</b>	Serum, Plasma, Tissue Homogenates, Cell Lysates, Cell Culture Supernates or Other Biological Fluids
<b>Uniprot ID</b>	<b>035622</b>
<b>Sensitivity</b>	12.4 pg/ml
<b>Entrez</b>	<b>14170</b>
<b>Expiration Date</b>	Please enquire.



Representative  
Standard  
Curve