


Product Datasheet

Thyroxine ELISA Kit (orb1473279)

Description	Thyroxine ELISA Kit
Reactivity	All
Range	1.25 ng/ml - 80 ng/ml
Tested Applications	EIA
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma or Other Biological Fluids
Sample Types	Serum, Plasma or Other Biological Fluids
Sensitivity	0.63 ng/ml
Expiration Date	Please enquire.



Representative
Standard
Curve