


Product Datasheet

Testosterone ELISA Kit (orb1473280)

Description	Testosterone ELISA Kit
Reactivity	All
Range	0.156 ng/ml - 10 ng/ml
Tested Applications	EIA
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Urine or Other Biological Fluids
Sample Types	Serum, Plasma, Urine or Other Biological Fluids
Sensitivity	0.063 ng/ml
Expiration Date	Please enquire.



Representative
Standard
Curve