


## Product Datasheet

### Estriol ELISA Kit (orb1473283)

<b>Description</b>	Estriol ELISA Kit
<b>Reactivity</b>	All
<b>Range</b>	1.56 ng/ml - 100 ng/ml
<b>Tested Applications</b>	EIA
<b>Storage</b>	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Tissue Homogenates, Cell Lysates, Cell Culture Supernates or Other Biological Fluids
<b>Sample Types</b>	Serum, Plasma, Tissue Homogenates, Cell Lysates, Cell Culture Supernates or Other Biological Fluids
<b>Sensitivity</b>	0.52 ng/ml
<b>Expiration Date</b>	Please enquire.



Representative  
Standard  
Curve