


Product Datasheet

Estradiol ELISA Kit (orb1473284)

Description	Estradiol ELISA Kit
Reactivity	All
Range	0.156 ng/ml - 10 ng/ml
Tested Applications	EIA
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Serum, Plasma, Tissue Homogenates, Cell Lysates, Cell Culture Supernates or Other Biological Fluids
Sample Types	Serum, Plasma, Tissue Homogenates, Cell Lysates, Cell Culture Supernates or Other Biological Fluids
Sensitivity	0.062 ng/ml
Expiration Date	Please enquire.



Representative Standard Curve