


Product Datasheet

Swine IL-8 ELISA Kit (orb1473285)

| | |
|----------------------------|---|
| Description | Swine IL-8 ELISA Kit |
| Reactivity | Porcine |
| Range | 62.5 pg/ml - 4000 pg/ml |
| Tested Applications | ELISA |
| Target | IL8 |
| Storage | Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Application notes | Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma |
| Sample Types | Cell Culture Supernates, Serum, Plasma |
| Uniprot ID | P26894 |
| Sensitivity | 30 pg/ml |
| Entrez | 396880 |
| Expiration Date | Please enquire. |



Representative Standard Curve