


## Product Datasheet

### Swine IL-1 alpha/IL-1F1 ELISA Kit (orb1473288)

<b>Description</b>	Swine IL-1 alpha/IL-1F1 ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	93.6 pg/ml - 6000 pg/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	IL1A
<b>Storage</b>	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
<b>Sample Types</b>	Cell Culture Supernates, Serum, Plasma
<b>Uniprot ID</b>	<b>P18430</b>
<b>Sensitivity</b>	50 pg/ml
<b>Entrez</b>	<b>397094</b>
<b>Expiration Date</b>	Please enquire.



Representative Standard Curve