


Product Datasheet

Rat Galectin 1 ELISA Kit (orb1473314)

Description	Rat Galectin 1 ELISA Kit
Reactivity	Rat
Range	156 pg/ml - 10000 pg/ml
Tested Applications	ELISA
Target	LGALS1
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P11762
Sensitivity	80 pg/ml
Entrez	56646
Expiration Date	Please enquire.



Representative Standard Curve