


Product Datasheet

Mouse PLA2G2A ELISA Kit (orb1473353)

Description	Mouse PLA2G2A ELISA Kit
Reactivity	Mouse
Range	31.2 pg/ml - 2000 pg/ml
Tested Applications	ELISA
Target	PLA2G2A
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	P31482
Sensitivity	15 pg/ml
Entrez	18780
Expiration Date	Please enquire.



Representative Standard Curve