


Product Datasheet

Mouse LI Cadherin ELISA Kit (orb1473360)

Description	Mouse LI Cadherin ELISA Kit
Reactivity	Mouse
Range	46.9 pg/ml - 3000 pg/ml
Tested Applications	ELISA
Target	CDH17
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9R100
Sensitivity	20 pg/ml
Entrez	12557
Expiration Date	Please enquire.



Representative
Standard
Curve