


Product Datasheet

Mouse EGFL7 ELISA Kit (orb1473375)

Description	Mouse EGFL7 ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4000 pg/ml
Tested Applications	ELISA
Target	EGFL7
Storage	Shipped and store at 4°C for 6 months, store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Application notes	Detection Method: Colorimetric, Sample Types: Cell Culture Supernates, Serum, Plasma
Sample Types	Cell Culture Supernates, Serum, Plasma
Uniprot ID	Q9QXT5
Sensitivity	30 pg/ml
Entrez	353156
Expiration Date	Please enquire.



Representative Standard Curve